

# Create Project from Pairwise Analysis

Any gene list generated in **Pairwise Analysis** can be saved and further analyzed using additional features in **Project Analysis**.

1. Select **Save as Project** to from **Pairwise Analysis**.
2. Provide a **Title** and **Description** for the new (sub)project.
3. Provide optional **Notes**, for example, how the gene list was generated.

Pairwise Analysis: mouse [ Scatter Plot | Export Results | Ontology | **1** | Save As Project ]

	Group 1	Group 2
Conditions:	Day 0	Day 2
Experiments:	23917, 23918	23919, 23920
Significance:	10, None	
Normalization:	Mean Intensity	
Quality Cutoff:	None	
Data Transformation:	None	

Show: 20 Search (38 results found) [1 - 20] [21 - 38]

GeneSifter.Net - Microarray Analysis Online - Microsoft Internet Explorer

Address: https://gs1.genesifter.net/users/?action=sub\_project\_info\_pw&mod=arr&file=253.txt&log=2&tag

GeneSifter.Net™ Home | About vizXlabs

**Control Panel**

- Analysis
  - Pairwise
  - Projects
- Inventories
  - Projects
  - Experiments
  - Arrays
  - Conditions
  - Targets
  - Flexload Protocols
  - MT&MF

**2**

**3**

**>> Sub Project Information**

Title:

Description:

Notes:

Save Sub Project Cancel

https://gs1.genesifter.net/users/ Internet